

BPO021-06

会場: 301B

時間: 5月25日10:15-10:30

MC-ICP-MSと気化法を用いたRe-Os同位体迅速測定法の開発：古海洋環境変動の解明に向けて

A rapid determination method for Re and Os isotope compositions using MC-ICP-MS combined with sparging method

野崎 達生^{1*}, 鈴木 勝彦¹, G. Ravizza², 木村 純一¹, 常青¹

Tatsuo Nozaki^{1*}, Katsuhiko Suzuki¹, Gregory Ravizza², Jun-Ichi Kimura¹, Qing Chang¹

¹(独)海洋研究開発機構, ²ハワイ大学

¹JAMSTEC/IFREE, ²Univ. of Hawaii

We present a rapid determination method for Re and Os isotope compositions using MC-ICP-MS (NEPTUNE) combined with a sparging method. The sparging method allows us to measure Os isotope compositions just after acid sample digestion in a Carius tube without the most commonly used but time-consuming solvent extraction by using CCl₄ and HBr. Our NEPTUNE is equipped with the multi-ion counter (MIC) system, enabling simultaneous measurement of a maximum of four Os isotopes. Compared to negative thermal ionization mass spectrometry (N-TIMS) which is now the most widely used for Re-Os analyses, sample throughput becomes several times higher. Therefore, the Re and Os measurement of sedimentary rock samples by NEPTUNE combined with the sparging method is expected to be a powerful tool to reconstruct the secular change of the marine Os isotope compositions with high resolution, which possibly unravels the cause of paleo-ocean global environmental change. In this presentation, we introduce the developed analytical method of Re and Os isotopes and their geological application. The accuracy, reproducibility and Re-Os data of GSJ geochemical reference samples (JCh-1 and JMS-2) measured by NEPTUNE are also compared with those by N-TIMS.

キーワード: Re-Os同位体, MC-ICP-MS, NEPTUNE, 気化法, マルチイオンカウンター(MIC)

Keywords: Re-Os isotope, MC-ICP-MS, NEPTUNE, sparging method, multi-ion counter (MIC)