

SEM001-P02

Room:Convention Hall

Time:May 26 10:30-13:00

Ultrafine-scale magnetostratigraphy with SQUID microscope: Application to ferromanganese crust and other materials

Hirokuni Oda^{1*}, Isoji MIYAGI¹, Akira Usui¹, Benjamin P. Weiss³, Franz J. Baudenbacher⁴, Eduardo A. Lima³

¹Geological Survey of Japan, AIST, ²Kochi University, ³Massachusetts Institute of Technology, ⁴Vanderbilt University

Scanning SQUID microscopy enables us to do mapping of magnetic field over thin sections of geological samples at sub-millimeter scale. With this new technique, ultrafine magnetostratigraphy can be achieved on geological sample which we could not observe the polarity pattern due to the limitation of the spacial resolution. In the presentation, we show a successful example on marine ferromanganese crust and extend the possibility in the future. Hydrogenetic ferromanganese crusts are iron-manganese oxide chemical precipitates on the seafloor that grow over periods of tens of millions of years. Their secular records of chemical, mineralogical, and textural variations are archives of deep-sea environmental changes. However, environmental reconstruction requires reliable high-resolution age dating. Earlier chronological methods using radiochemical and stable isotopes provided age models for ferromanganese crusts, but have limitations on the millimeter scale. For example, the reliability of $^{10}\text{Be}/^{9}\text{Be}$ chronometry, commonly considered the most reliable technique, depends on the assumption that the production and preservation of ^{10}Be are constant, and requires accurate knowledge of the ^{10}Be half-life. To overcome these limitations, we applied an alternative chronometric technique, magnetostratigraphy, to a 50-mm-thick hydrogenetic ferromanganese crust (D96-m4) from the northwest Pacific. Submillimeter-scale magnetic stripes originating from approximately oppositely magnetized regions oriented parallel to bedding were clearly recognized on thin sections of the crust using a high-resolution magnetometry technique called scanning SQUID (superconducting quantum interference device) microscopy. By correlating the boundaries of the magnetic stripes with known geomagnetic reversals, we determined an average growth rate of 5.1 ± 0.2 mm/m.y., which is within 16% of that deduced from the $^{10}\text{Be}/^{9}\text{Be}$ method (6.0 ± 0.2 mm/m.y.). This is the finest-scale magnetostratigraphic study of a geologic sample to date. Ultrafine-scale magnetostratigraphy using SQUID microscopy is a powerful new chronological tool for estimating ages and growth rates for hydrogenetic ferromanganese crusts. It provides chronological constraints with the accuracy promised by the astronomically calibrated magnetostratigraphic time scale (1-40 k.y.). The technique can be extended to other geological objects such as stalagmite, hydrothermal deposits, desert varnish, etc. An analogue of desert varnish would be found on Mars, which might have recorded ancient Martian magnetic field.